

INNOVACIÓN

N 1108 PROGRESSIVE

Apply different therapies and programs on 2 MUSCLE GROUPS at the SAME TIME...



RUSSIANS + SQUARE + T.E.N.S. + INTERFERENTIAL



EQUIPO BIFÁSICO · CORPORAL
54
PROGRAMAS

2 YEARS WARRANTY
+ TRAINING
+ KRYALING

Professional Equipment

- BIPHASIC CURRENTS.
- 8 INDEPENDENT CHANNELS.
- 2 SIMULTANEOUS WAVES.
- 2 DIGITAL SCREENS.
- PULSE WIDTH CONTROL.
- EMISSION STARTS and ENDS ON RAMP.
- 100% CONDUCTIVE ELECTRODES.
- 34 PRESET PROGRAMS.
- 20 CUSTOM PROGRAMS.

N 1108 PROGRESSIVE - DOUBLE THERAPY

N 1108 PROGRESSIVE® is an advanced electrotherapy device for esthetics, physiotherapy, sports, and rehab.

It's versatile, work with its eight output channels or select to perform sessions with four channels on one muscle group and the remaining four on other muscle group - simultaneously - and with different waves and programs.

This feature (unique in the market) allows professionals to work in different areas of the body at the same time, focusing on each muscle need and required therapy.

SMT MICROELECTRONICS

Its eight channels operate independently or together, each one with an isolated transformer.

N 1108 PROGRESSIVE® grants the possibility of delivering visible, effective, and durable results by utilizing Russians, Square, or T.E.N.S. waves, with **start and end** of its pulse emission **on ramp** and along its 54 programs (specially developed by professionals).

Modify all parameters in a **safe and rapid way**. Also, trace each patient evolution by creating and saving new programs in its memories.

NOTO
TRUSTED TECHNOLOGY

NOTO proposes a different course for aesthetics, a new scenario where therapies combination, cutting-edge technology, and quality products and supplies are fundamental for visible, effective, and durable results.

RUSSIAN WAVES

Are modulated medium frequency currents, delivered in low frequency signals. Its medium frequency objective is to decrease skin's impedance to allow low frequency work on muscles. These currents main target is muscle tone or body outline contouring. N 1108 is available to work with bipolar or tetrapolar interferential wave mode.

MAIN TREATMENTS:

- Flaccidity and body contouring.
- Glutei hypertrophy and tone.
- Progressive muscle hypertrophy.
- Atrophy by inactivity.
- Post surgery strengthen process.
- Sports rehab.

T.E.N.S.

T.E.N.S. programs utilizes nerve system mechanisms to mitigate chronic and acute pain.

Decrease of pain, creating an immediate relief without side effects, is T.E.N.S. main advantage.

MAIN TREATMENTS:

- Peripheral nerves injuries.
- Nerve compression and reflex sympathetic dystrophy (carpal tunnel syndrome).
- Stump and/or phantom member pain.
- Posttherpetic neuralgia.
- Back and neck pain.

SQUARE WAVES

They are presented in low and medium frequency pulses. These waves are the most energetic waves of all electrostimulation currents.

MAIN TREATMENTS:

- Sports.
- Body Contouring.
- Tone.
- Strengthening.
- Strength & endurance.
- Muscle firmnes.
- Contouring.

The device N 1108 comes with:



8 LARGE Electrodes

8 SMALL Electrodes

6 Elastic Straps

8 Connection cables


1 User's manual

Technical Data N 1108:

| | |
|-----------------------|-------------------------------------|
| MANUFACTURED BY: | NOTO GROUP S.A. |
| COUNTRY: | ARGENTINA |
| NAME OF PRODUCT: | N 1108 PROGRESSIVE |
| WIDE: | 350 MM |
| LENGHT: | 310 MM |
| HEIGHT: | 160 MM-75 MM |
| WEIGHT: | 3,1 KG |
| POWER SUPPLY VOLTAGE: | 98 to 240 VAC/60 HZ |
| POWER: | 21 WATT |
| PROGRAMMING: | MICROCONTROLLER |
| TYPES OF WAVES: | RUSSIANS, SQUARE, AND T.E.N.S. |
| MAXIMUM CURRENT: | 120 MA. Impedance 1K |
| QUANTITY OF CHANNELS: | 8 for Body, ISOLATED INDEPENDENTLY. |

| | |
|---------------|------------------|
| PROGRAMS: | 34 PREDETERMINED |
| TIMER: | 30 MIN. MAX. |
| CONTRACTION: | 2 to 15 SECONDS |
| RELAXATION: | 2 to 15 SECONDS |
| FREQUENCY: | 2 to 150 HZ |
| RELATION: | 1:1 - 2:1 - 4:1 |
| ASCENT RAMP: | 1 SECOND |
| DESCENT RAMP: | 1 SECOND |
| SEQUENCE: | 1,2,3,4 or 5 |

NOTE: DEVICE SPECIFICATIONS MAY CHANGE WITHOUT PREVIOUS NOTICE.



2500 HZ. RUSSIAN WAVES FREQUENCY CARRIER

NOTO supports all clients with a set of tools available at all times.

Tools as training, post-sale service support, and assured spare parts. All these instruments with the aim of making Noto's products your best choice.

MORE NOTO EQUIPMENTS:

N 1800 - MICRODERMABRASION + ELECTROPORATION WITH "ID" SYSTEM.
BIMAX - BIPOLAR RADIOFREQUENCY WITH CROSS FIELDS EMISSION.